RFID for medical, healthcare and analytics

Wireless ID and sensors

for item management, tracing and tracking, probe management, condition and quality monitoring

microsensys - your partner for wireless identification and sensors



Your Applications



Our Solution for



Surgical

- High temperature and pressure tolerance
- Viable through several sterilization processes, high temperature tolerance
- Miniaturized design for unobstructed use

Instrument

- On-metal (and non-metal) application or integration
- Internal memory for ID, life cycle management, documentation of use(s)
- Brand protection and certificate of authenticity

Tagging

- Compatible stationary or mobile readers
- ❖ Data transfer via HF (NFC), UHF
- Globally unique serial number for exact allocation
- Keep track via custom software or integration into existing system

Products:

Q2.5 or Q1.6 mini
TAG
iID® mobile readers
TELID®soft
laser printing option



Our Solution for



Sterilization

- Rugged design for harsh environments, high temperature and pressure tolerance
- ❖ Passive components with near unlimited life span

Container

- ❖ Simple and quick installation of TAGs and/ or labels
- On-metal (and non-metal) application or integration
- up to 16 kbit RW memory for ID or description

Identification

- Compatible stationary or mobile readers
- ❖ Data transfer via HF or UHF
- ❖ Globally unique serial number for exact allocation
- Cataloguing by custom software or integration into existing system

Products:

Q2.5 mini TAG
DUALTAGS
On-metal Plates
iID® mobile readers
TELID®soft
Barcode or matrix
code laser printing
option

Allocation and Cataloguing



Our Solution for



Wireless

- Simple and quick installation in hard-to-reach areas
- On-metal (and non-metal) application or integration
- Data transfer via HF, UHF or Bluetooth Low Energy

Process

- Continuous and/or real-time supervision of temperature, pressure, humidity and more
- Rugged design for harsh environments, high temperature and pressure tolerance

Monitoring

- Freely adjustable parameters and limits
- Compatible stationary or mobile readers
- Analysis by custom software or integration into existing system

Products:

TELID®3T
iID® Desktop
iID® mobile readers
TELID®station
TELID®soft

General Operation

Process circle from initialization to evaluation



INITIALIZE

- Adjust sensor parameters
- •Adjust LIMIT CONTROL



Your process

EVALUATE

- •Read entire data
- •Evaluate and create graphical REPORTs
- •Save and export datasets for extended evaluation

PROGRAM

- •Set individual process **PARAMETERs**
- Add optional REMARKs
- Start data logging



10


CHECK

- •Use unique TELID® ID NUMBER
- •Check process quality using LIMIT CONTROL
- •Supported by smart mobile devices



System Overview Stationary - RFID sensor data loggers


TELID®3T solution comes with different kind of data loggers, read/write devices and TELID®soft or TELID®service developed for different operation systems. TELID®bundle unites all the components in one.



- 1) Cloud based data exchange on request in OEM projects
- 2) TELID®library is available and allows TELID® system integration into own infrastructure

System Overview Mobile - NFC sensor data loggers

TELID®3T NFC is built on TELID®312.nfc data logger & TELID®soft and can be used with microsensys read/write devices or **android smart phones** with NFC communication.



¹⁾ Cloud based data exchange on request in OEM projects

²⁾ TELID®library is available for Android and Windows and allows TELID® system integration into own infrastructure



TELID® Bundle Options

- ❖ TELID® Bundles: receive a complete bundle of several sensor data loggers, a compatible mobile or stationary reader, all necessary accessories and the base configuration of our software
- Plug-and-play system solution
- Rental possible for testing purposes

Not quite what you're looking for?

All our bundles can be customized and tailored to your needs. Still unsure? Don't worry! Our equipment can be rented for testing purposes and adjusted to your application. Please contact our sales team for more information our check out our website: www.microsensys.de

Technical Information - TELID®311, 314, 343

TELID [®] Type Product Code	Picture	Application Field	Compatible Readers	Technical
311 14.311.484.06	TELID® 311 VID conditions buyer 30°C \$125°C U battery insides (1 battery insides	Disinfecting, washing and drying processes +100°C	ilD®DESKTOPsmart ilD®POCKETwork ilD®M30 HEAD	dimensions: D27mm, thickness 12mm environmental condition: working temp. up to +100°C samples: up to 8000 sample time: 10s 256min
343.HTK 14.343.439.103		Disinfecting and washing processes +100°C	ilD®DESKTOPsmart ilD®POCKETwork ilD®M30 HEAD	dimensions: G¼" fitting, SW22, max length 48mm pressure range: 0 Bar 3Bar / 100Bar environmental condition: working temp. up to +100°C samples: up to 8000 sample time: 1min 59min
311.ac 14.311.489.00	TELIO® 311.8C artic temperature topor - SCC - +480°C (Li butterly finder)	Steam sterilization autoclave +135°C	iID®DESKTOPsmart iID®POCKETwork iID®M30 HEAD	dimensions: D27mm, thickness 12mm environmental condition: working temp. up to +140°C absolute pressure up to 5bar samples: up to 8000 sample time: 1s 59min
314.ac 14.314.499.01 14.314.499.05 14.314.499.07		Steam sterilization autoclave +135°C	iID®DESKTOPsmart iID®POCKETwork iID®M30 HEAD	dimensions: D27mm, thickness 12mm external sensor: pin length 15, 50 and 100mm environmental condition: working temp. up to +140°C absolute pressure up to 5bar samples: up to 8000 sample time: 1s 59min



Micro-Sensys GmbH

In der Hochstedter Ecke 2 D 99098 Erfurt Germany

TEL +49 361 59874 0 FAX +49 361 59874 17 WEB <u>www.microsensys.de</u>

microsensys US - for customer specific RFID solutions

Office Gainesville 4111 NW 6th Street, Suit D Gainesville, FL 32609 Florida, USA

Sales Office

150 North Michigan Avenue, 35th Floor Chicago, III 60601 Illinois, USA

WEB www.microsensys.us